

## Product datasheet for **GA109030**

### Human Zonula occludens protein 3 (TJP3) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	TJP3 CRISPRa kit - CRISPR gene activation of human tight junction protein 3
Symbol:	TJP3
UniProt ID:	<a href="#">O95049</a>
Synonyms:	ZO-3; ZO3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)