

## Product datasheet for **GA109040**

### Human ZNF777 activation kit by CRISPRa

#### Product data:

|                      |                                                                                                                                                                                                                                                                                                               |
|----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type:</b> | CRISPRa Kits                                                                                                                                                                                                                                                                                                  |
| <b>Description:</b>  | ZNF777 CRISPRa kit – CRISPR gene activation of human zinc finger protein 777                                                                                                                                                                                                                                  |
| <b>Symbol:</b>       | ZNF777                                                                                                                                                                                                                                                                                                        |
| <b>Components:</b>   | GA109040G1, ZNF777 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA109040G2, ZNF777 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA109040G3, ZNF777 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| <b>UniProt ID:</b>   | <u><a href="#">Q9ULD5</a></u>                                                                                                                                                                                                                                                                                 |
| <b>Synonyms:</b>     | KIAA1285                                                                                                                                                                                                                                                                                                      |
| <b>Format:</b>       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                       |

#### Product images:

