

Product datasheet for **GA109050**

Human **SLC39A1** activation kit by **CRISPRa**

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | SLC39A1 CRISPRa kit - CRISPR gene activation of human solute carrier family 39 member 1 |
| Symbol: | SLC39A1 |
| UniProt ID: | Q9NY26 |
| Synonyms: | ZIP1; ZIRTL |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)