

Product datasheet for **GA109053**

Human IL37 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	IL37 CRISPRa kit - CRISPR gene activation of human interleukin 37
Symbol:	IL37
UniProt ID:	<u>Q9NZH6</u>
Synonyms:	FIL1; FIL1(ZETA); FIL1Z; IL-1F7; IL-1H; IL-1H4; IL-1RP1; IL-37; IL1F7; IL1H4; IL1RP1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)