

Product datasheet for **GA109054**

Human IL36 alpha (IL36A) activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | IL36A CRISPRa kit - CRISPR gene activation of human interleukin 36 alpha |
| Symbol: | IL36A |
| UniProt ID: | Q9UHAZ |
| Synonyms: | FIL1; FIL1(EPSILON); FIL1E; IL-1F6; IL1(EPSILON); IL1F6 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)