

Product datasheet for **GA109076**

Human SIT (SIT1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SIT1 CRISPRa kit - CRISPR gene activation of human signaling threshold regulating transmembrane adaptor 1
Symbol:	SIT1
UniProt ID:	Q9Y3P8
Synonyms:	SIT; SIT-R
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)