

Product datasheet for **GA109088**

Human PCDH17 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | PCDH17 CRISPRa kit - CRISPR gene activation of human protocadherin 17 |
| Symbol: | PCDH17 |
| UniProt ID: | Q14917 |
| Synonyms: | PCDH68; PCH68 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)