

## Product datasheet for **GA109092**

### Human LSM3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	LSM3 CRISPRa kit - CRISPR gene activation of human LSM3 homolog, U6 small nuclear RNA and mRNA degradation associated
Symbol:	LSM3
UniProt ID:	<a href="#">P62310</a>
Synonyms:	SMX4; USS2; YLR438C
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)