

Product datasheet for **GA109093**

Human TINAG activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TINAG CRISPRa kit - CRISPR gene activation of human tubulointerstitial nephritis antigen
Symbol:	TINAG
UniProt ID:	<u>Q9UJW2</u>
Synonyms:	TIN-AG
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)