


Product datasheet for **GA109101**

Human GIDRP88 (R3HCC1L) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	R3HCC1L CRISPRa kit - CRISPR gene activation of human R3H domain and coiled-coil containing 1 like
Symbol:	R3HCC1L
Components:	GA109101G1, GIDRP88 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA109101G2, GIDRP88 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA109101G3, GIDRP88 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q7Z5L2
Synonyms:	C10orf28; GIDRP86; GIDRP88; PSORT
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)