

Product datasheet for **GA109103**

Human SMPDL3B activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SMPDL3B CRISPRa kit - CRISPR gene activation of human sphingomyelin phosphodiesterase acid like 3B
Symbol:	SMPDL3B
UniProt ID:	Q92485
Synonyms:	ASML3B
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)