

Product datasheet for **GA109106**

Human TP53TG5 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TP53TG5 CRISPRa kit - CRISPR gene activation of human TP53 target 5
Symbol:	TP53TG5
Components:	GA109106G1, TP53TG5 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA109106G2, TP53TG5 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA109106G3, TP53TG5 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9Y2B4</u>
Synonyms:	C20orf10; CLG01
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)