

Product datasheet for **GA109118**

Human hnRNP G (RBMX) activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | RBMX CRISPRa kit - CRISPR gene activation of human RNA binding motif protein X-linked |
| Symbol: | RBMX |
| UniProt ID: | P38159 |
| Synonyms: | hnRNP-G; HNRNPG; HNRPG; MRXS11; RBMXP1; RBMXRT; RNMX |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)