

Product datasheet for **GA109121**

Human **GW182 (TNRC6A)** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	TNRC6A CRISPRa kit - CRISPR gene activation of human trinucleotide repeat containing adaptor 6A
Symbol:	TNRC6A
UniProt ID:	Q8NDV7
Synonyms:	CAGH26; GW1; GW182; TNRC6
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)