

Product datasheet for **GA109128**

Human EIF3K activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	EIF3K CRISPRa kit - CRISPR gene activation of human eukaryotic translation initiation factor 3 subunit K
Symbol:	EIF3K
UniProt ID:	<u>Q9UBQ5</u>
Synonyms:	ARG134; EIF3-p28; EIF3S12; HSPC029; M9; MSTP001; PLAC-24; PLAC24; PRO1474; PTD001
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)