

Product datasheet for **GA109136**

Human ProSAAS (PCSK1N) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PCSK1N CRISPRa kit - CRISPR gene activation of human proprotein convertase subtilisin/kexin type 1 inhibitor
Symbol:	PCSK1N
Components:	GA109136G1, ProSAAS gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA109136G2, ProSAAS gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA109136G3, ProSAAS gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9UHG2</u>
Synonyms:	BigLEN; PEN; PROSAAS; SAAS; SCG8; SgVIII
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)