

Product datasheet for **GA109137**

Human **KCNMB4** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	KCNMB4 CRISPRa kit - CRISPR gene activation of human potassium calcium-activated channel subfamily M regulatory beta subunit 4
Symbol:	KCNMB4
UniProt ID:	Q86W47
Synonyms:	calcium-activated potassium channel beta 4 subunit; large conductance calcium-dependent potassium ion channel beta 4 subunit; potassium large conductance calcium-activated channel, subfamily M, beta member 4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

