

## Product datasheet for **GA109154**

### Human PPP2R3B activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	PPP2R3B CRISPRa kit - CRISPR gene activation of human protein phosphatase 2 regulatory subunit B"beta
Symbol:	PPP2R3B
UniProt ID:	<a href="#">Q9Y5P8</a>
Synonyms:	NYREN8; PPP2R3L; PPP2R3LY; PR48
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)