

## Product datasheet for **GA109167**

### Human robld3 (LAMTOR2) activation kit by CRISPRa

#### Product data:

|               |   |
|---------------|---|
| Product Type: | CRISPRa Kits  |
| Description:  | LAMTOR2 CRISPRa kit - CRISPR gene activation of human late endosomal/lysosomal adaptor, MAPK and MTOR activator 2 |
| Symbol:       | LAMTOR2   |
| UniProt ID:   | <a href="#">Q9Y2Q5</a>  |
| Synonyms:     | ENDAP; HSPC003; MAPBPIP; MAPKSP1AP; p14; Ragulator2; ROBLD3   |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)   |

#### Product images:



[View online >](#)