

Product datasheet for **GA109168**

Human MRPS28 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	MRPS28 CRISPRa kit - CRISPR gene activation of human mitochondrial ribosomal protein S28
Symbol:	MRPS28
Components:	GA109168G1, MRPS28 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA109168G2, MRPS28 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA109168G3, MRPS28 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9Y2Q9</u>
Synonyms:	HSPC007; MRP-S28; MRP-S35; MRPS35
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)