

Product datasheet for **GA109201**

Human **SLC43A3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	SLC43A3 CRISPRa kit - CRISPR gene activation of human solute carrier family 43 member 3
Symbol:	SLC43A3
UniProt ID:	<u>Q8NBI5</u>
Synonyms:	EEG1; FOAP-13; PRO1659; SEEEG-1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)