

## Product datasheet for **GA109205**

### Human C20orf30 (TMEM230) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	TMEM230 CRISPRa kit - CRISPR gene activation of human transmembrane protein 230
Symbol:	TMEM230
UniProt ID:	<a href="#">Q96A57</a>
Synonyms:	C20orf30; dj1116H23.2.1; HSPC274
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)