

Product datasheet for **GA109206**

Human WDR91 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	WDR91 CRISPRa kit - CRISPR gene activation of human WD repeat domain 91
Symbol:	WDR91
Components:	GA109206G1, WDR91 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA109206G2, WDR91 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA109206G3, WDR91 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>A4D1P6</u>
Synonyms:	HSPC049; SORF-1; SORF1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)