

Product datasheet for **GA109210**

Human **CCDC113** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	CCDC113 CRISPRa kit - CRISPR gene activation of human coiled-coil domain containing 113
Symbol:	CCDC113
Components:	GA109210G1, CCDC113 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA109210G2, CCDC113 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA109210G3, CCDC113 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9H0I3</u>
Synonyms:	HSPC065
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)