


## Product datasheet for **GA109214**

### Human NDUFAF4 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                                  |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                     |
| Description:  | NDUFAF4 CRISPRa kit - CRISPR gene activation of human NADH:ubiquinone oxidoreductase complex assembly factor 4                                                                                                                                                                                                   |
| Symbol:       | NDUFAF4                                                                                                                                                                                                                                                                                                          |
| Components:   | GA109214G1, NDUFAF4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA109214G2, NDUFAF4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA109214G3, NDUFAF4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <a href="#">Q9P032</a>                                                                                                                                                                                                                                                                                           |
| Synonyms:     | bA22L21.1; C6orf66; HRPAP20; HSPC125; My013                                                                                                                                                                                                                                                                      |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                          |

#### Product images:



[View online »](#)