

Product datasheet for **GA109239**

Human Secretogranin 3 (SCG3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SCG3 CRISPRa kit - CRISPR gene activation of human secretogranin III
Symbol:	SCG3
UniProt ID:	<u>Q8WXD2</u>
Synonyms:	SGIII
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)