

## Product datasheet for **GA109245**

### Human TAGLN3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	TAGLN3 CRISPRa kit - CRISPR gene activation of human transgelin 3
Symbol:	TAGLN3
UniProt ID:	<a href="#">Q9UI15</a>
Synonyms:	NP22; NP24; NP25
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)