

Product datasheet for **GA109261**

Human PACSIN3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PACSIN3 CRISPRa kit - CRISPR gene activation of human protein kinase C and casein kinase substrate in neurons 3
Symbol:	PACSIN3
UniProt ID:	Q9UKS6
Synonyms:	SDPIII
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)