

Product datasheet for **GA109263**

Human Tropomodulin 3 (TMOD3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMOD3 CRISPRa kit – CRISPR gene activation of human tropomodulin 3
Symbol:	TMOD3
Components:	GA109263G1, Tropomodulin 3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA109263G2, Tropomodulin 3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA109263G3, Tropomodulin 3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9NYL9</u>
Synonyms:	UTMOD
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

