

Product datasheet for **GA109294**

Human PSMC3IP activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PSMC3IP CRISPRa kit - CRISPR gene activation of human PSMC3 interacting protein
Symbol:	PSMC3IP
UniProt ID:	<u>Q9P2W1</u>
Synonyms:	GT198; HOP2; HUMGT198A; ODG3; TBPIP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)