

Product datasheet for **GA109299**

Human **SAC3D1** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	SAC3D1 CRISPRa kit - CRISPR gene activation of human SAC3 domain containing 1
Symbol:	SAC3D1
Components:	GA109299G1, SAC3D1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA109299G2, SAC3D1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA109299G3, SAC3D1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	HSU79266; SHD1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)