

Product datasheet for **GA109310**

Human RMC1 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | RMC1 CRISPRa kit - CRISPR gene activation of human regulator of MON1-CCZ1 |
| Symbol: | RMC1 |
| UniProt ID: | Q96DM3 |
| Synonyms: | HsT2591; Mic-1; MIC1 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)