

## Product datasheet for **GA109326**

### Human PKN beta (PKN3) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	PKN3 CRISPRa kit - CRISPR gene activation of human protein kinase N3
Symbol:	PKN3
UniProt ID:	<a href="#">Q6P5Z2</a>
Synonyms:	UTDP4-1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)