


Product datasheet for **GA109327**

Human PURG activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PURG CRISPRa kit - CRISPR gene activation of human purine rich element binding protein G
Symbol:	PURG
Components:	GA109327G1, PURG gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA109327G2, PURG gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA109327G3, PURG gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9UJV8
Synonyms:	PURG-A; PURG-B
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)