

## Product datasheet for **GA109329**

### Human PNMA3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	PNMA3 CRISPRa kit - CRISPR gene activation of human PNMA family member 3
Symbol:	PNMA3
UniProt ID:	<a href="#">Q9UL41</a>
Synonyms:	MA3; MA5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)