

Product datasheet for GA109331

Human SERTAD3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SERTAD3 CRISPRa kit – CRISPR gene activation of human SERTA domain containing 3
Symbol:	SERTAD3
Components:	GA109331G1, SERTAD3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GEI00074) GA109331G2, SERTAD3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GEI00074) GA109331G3, SERTAD3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GEI00074) 1 CRISPRa-Enhancer vector, SKU GEI00056 1 CRISPRa scramble vector, SKU GEI00077
UniProt ID:	<u>Q9UJW9</u>
Synonyms:	RBT1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

