

Product datasheet for **GA109336**

Human PDZRN4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PDZRN4 CRISPRa kit - CRISPR gene activation of human PDZ domain containing ring finger 4
Symbol:	PDZRN4
UniProt ID:	<u>Q6ZMN7</u>
Synonyms:	LN4; SAMCAP3L
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)