

Product datasheet for **GA109365**

Human PACSIN1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PACSIN1 CRISPRa kit - CRISPR gene activation of human protein kinase C and casein kinase substrate in neurons 1
Symbol:	PACSIN1
UniProt ID:	Q9BY11
Synonyms:	SDPI
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

