

Product datasheet for **GA109376**

Human SH3KBP1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SH3KBP1 CRISPRa kit - CRISPR gene activation of human SH3 domain containing kinase binding protein 1
Symbol:	SH3KBP1
UniProt ID:	Q96B97
Synonyms:	CD2BP3; CIN85; GIG10; HSB-1; HSB1; MIG18
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)