

Product datasheet for **GA109388**

Human CSEN (KCNIP3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KCNIP3 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel interacting protein 3
Symbol:	KCNIP3
UniProt ID:	Q9Y2W7
Synonyms:	CSEN; DREAM; KCHIP3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)