

Product datasheet for **GA109389**

Human KCHIP2 (KCNIP2) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KCNIP2 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel interacting protein 2
Symbol:	KCNIP2
UniProt ID:	<u>Q9NS61</u>
Synonyms:	KCHIP2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)