

Product datasheet for **GA109390**

Human KCNIP1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KCNIP1 CRISPRa kit - CRISPR gene activation of human potassium voltage-gated channel interacting protein 1
Symbol:	KCNIP1
UniProt ID:	Q9NZI2
Synonyms:	KCHIP1; VABP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)