

Product datasheet for **GA109392**

Human ZNF354C activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF354C CRISPRa kit - CRISPR gene activation of human zinc finger protein 354C
Symbol:	ZNF354C
UniProt ID:	<u>Q86Y25</u>
Synonyms:	KID3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)