

Product datasheet for **GA109402**

Human PI 3 Kinase regulatory subunit 4 (PIK3R4) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PIK3R4 CRISPRa kit - CRISPR gene activation of human phosphoinositide-3-kinase regulatory subunit 4
Symbol:	PIK3R4
UniProt ID:	<u>Q99570</u>
Synonyms:	p150; VPS15
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)