

Product datasheet for **GA109410**

Human WDR8 (WRAP73) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	WRAP73 CRISPRa kit - CRISPR gene activation of human WD repeat containing, antisense to TP73
Symbol:	WRAP73
UniProt ID:	Q9P2S5
Synonyms:	WDR8
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

□



[View online »](#)