

Product datasheet for **GA109427**

Human UBQLN3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	UBQLN3 CRISPRa kit - CRISPR gene activation of human ubiquilin 3
Symbol:	UBQLN3
Components:	GA109427G1, UBQLN3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA109427G2, UBQLN3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA109427G3, UBQLN3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9H347</u>
Synonyms:	TUP-1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

