


Product datasheet for **GA109445**

Human WDR42A (DCAF8) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DCAF8 CRISPRa kit - CRISPR gene activation of human DDB1 and CUL4 associated factor 8
Symbol:	DCAF8
Components:	GA109445G1, WDR42A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA109445G2, WDR42A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA109445G3, WDR42A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5TAQ9</u>
Synonyms:	GAN2; H326; WDR42A
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)