

Product datasheet for **GA109446**

Human **KCNK4** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	KCNK4 CRISPRa kit - CRISPR gene activation of human potassium two pore domain channel subfamily K member 4
Symbol:	KCNK4
UniProt ID:	Q9NYG8
Synonyms:	K2p4.1; K2P4.1 potassium channel; potassium channel, subfamily K, member 4; TRAAK; TRAAK, DKFZP566E164; TRAAK1; TWIK-related arachidonic acid-stimulated potassium channel protein; two pore K ⁺ channel KT4.1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)