

Product datasheet for **GA109456**

Human **TAS2R3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	TAS2R3 CRISPRa kit - CRISPR gene activation of human taste 2 receptor member 3
Symbol:	TAS2R3
UniProt ID:	Q9NYW6
Synonyms:	T2R3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)