

Product datasheet for **GA109457**

Human **TAS2R4** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	TAS2R4 CRISPRa kit - CRISPR gene activation of human taste 2 receptor member 4
Symbol:	TAS2R4
UniProt ID:	<u>Q9NYW5</u>
Synonyms:	T2R4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)